

Product display 智能测肤

Inspired by:

The designer originally designed the Capricorn imager from the sci-fi movie aircraft. In addition, the designer himself has always yearned for a sense of science and technology, mystery, and exploring the unknown, so he openedCapricorn. The road to exploration of smart imagers.

Product details



28 MILLION HD PIXELS

Comparable to the camera of high-speed snapping instrument, with dear imaging, combined with 8 spectrum to quickly identify various skin problems



DISPLAY

2k high definition display, resolution up to 1920x1080, can be moved freely

COMBINED FAN LEAF SHADING SHEET

Effectively block the external light source, improve the image quality, reduce the analysis error.



3.0 HIGH-SPEED INFORMATION TRANSMISSION LINE

Blue high- speed transmission head, purecopper transmission line, to ensure the quality and speed of transmission information






SWITCH

Using SCR power switch, effective Stable voltage to protect the instrument

QANTITATIVE AND DIFFERENTIATED DATA OF 14 INDICATORS

www.sincoherengroup.com



SKIN DIAGNOSIS


The data of each indicator is presented separately by a professional dermatologist who can broadcast the quantitative skin measurement results by double click

Software Introduction


3D simulation technology: 3D simulation accurately recognizes facial features and quantitatively analyzes four major problem skins

Deep learning technology: Combine with large database annotations to accurately identify and analyze problems existing in various skins.

Ai face recognition technology: using artificial intelligence face recognition algorithm to accurately identify o area t area, Data is more accurate and faster.



ai intelligent skin detection 8 spectrum interpretation of different theme interfaces



Ten modules of Capri AI intelligent imager

the file management module Cloud computing and cloud storage modules Third, the surface and deep quantitative analysis module

Fourth, the skin prediction module Five, problem skin

module Sixth, 3d simulation module

Seven, comparative analysis module Eight,


comprehensive report module Nine,

solution module Ten, VIP custom

modules




www.sincoherengroup.com




Newly upgraded 8 spectra for surface, deep, and quantitative analysis of 14 major skin indicators






www.sincoherengroup.com

Four major problem skins combined with the total data of 14 indicators to quantify the differentiation






Three major indicators for predicting skin status, future aging, and skin beauty prediction



Combining total data to obtain indicators of comprehensive reports All pages can be scanned and printed for analysis reports





APPLICATION



Before and After



Before



After



Before



After



Before After



Before



After



Skin Analyzer *

SPECIFICATION

Illumination (RGB)	8,800 lux
Fluorescent tube life	9000 h
Color temperature	7200 K
Video output	NT
CPU/ Pentium 4	1.8G H z
Computer & Monitor	OS: Windows XP or Windows 7
N.W/G.W	13kg/15kg
Dimension	49cm*52cm*45cm (L×W×H)
Max Resolutions	10megapixel
Hard Disk Space	120GB
RAM	2GB
Power Device Electrical motor	48.5*44*50 cm
Camera	1:1.7"CCD Digital Camera
USB	2.0 Port
Electrical Requirement	AC 220V±10% 50Hz
Monitor 1280X800 pixels resolution	
with Power cart AC in 220v (or 110v	Switch able)
Environment Temperature	10-35 °C



Add: A-4, Sinotrans Plaza, No. 43 Xizhimen North St,Beijing 100044 C

Tel: +86 18701341926

Web: www.sincoherengroup.com







