





A beauty machine with multifuntion, non-invasive, lifting and contouring treatment

Illustrate of machine

1) Overall appearance

Side of the machine



Front of the machine



Illustrate of machine

2) HIFU hand-piece and Cartridges

1.5mm - 16.0mm







Illustrate of machine

3) Ultra system hand-piece and Cartridges



Vmax Handpiece

- ① Cartridge ultra 7-1.5mm
- ② Cartridge ultra4-3.0mm
- ③Cartridge ultra4-4.5mm
- 4 Cartridge ultra4-8.0mm
- ©Cartridge ultra4-13.0mm




Illustrate of machine

4) User Interface









Advantages



A. Features of the upgrade HIFU:

- 1) Quick & short treatment time: 30 MINS one face treatment
- 2) SMAS contraction : collagen remodeling , elastine fiber contraction
- 3) None downtime: skin just become red within the first several hours, then skin recover.
- 4) Instant result will be checked from the second month to the ninth month, good result will last 2-3 years
- 5) Totally non-invasive

B. Treatment advantages

- 1) Radar carving:1-9 points that are adjusting can be shot per second.
- 2) New radar technology cartridges are more suitable for some small areas, like areas around eyes, mouth and forehead.

Applications



- Removing wrinkles on around forehead, eyes, mouth, etc
- 2) Lifting and tightening both cheeks skin.
- 3) Improving skin elasticity and shaping contour.
- 4) Improving jaw line, reducing "marionette lines".
- 5) Tightening the skin tissue on forehead, lifting the eyebrows lines.
- 6) Improving skin complexion, making the skin delicate and bright.
- 7) Removing neck wrinkles, protecting neck aging.
- 8) Weight loss.

BEFORE & AFTER





After



Refore



After



Before



After



Before



After

Technical Parameter

Screen	15" color touch LCD screen
Cartridge frequency	4MHZ, 7MHZ
HIFU Energy	0.2-1.2J (0.1J step, 18 steps)
HIFU Length	5-25mm (1.0mm step, 20steps)
HIFU Equipped cartridge	1.5mm / 3.0mm / 4.5mm
HIFU Optional cartridge	8.0mm / 6mm / 10mm / 13mm
	/16mm
Vmax Equipped cartridge	1.5mm / 3.0mm / 4.5mm
Vmax Optional cartridge	Breast 4.5mm / 8.0mm
	/13.0mm
Lifespan of Probe	62000 shots/Probe(Vmax)
	20000 shots/probe(HIFU)
Package Size	54cm*55cm*45cm
Gross weight	12kg
Voltage	AC110V-240V.50/60Hz

